- DE Gebrauchsanweisung Nephroskope
- EN Instruction for Use Nephroscopes
- ES Instrucciones de uso de los nefroscopios
- FR Mode d'emploi des néphroscopes
- IT Istruzioni per l'uso dei nefroscopi

Bo87

- DE Diese Gebrauchsanweisung betrifft:
- EN These instructions for use concern:
- ES Estas instrucciones de uso se refieren a:
- FR Ce mode d'emploi concerne:
- IT Queste istruzioni per l'uso riguardano:

E.9600.00 E.9601.00 E.9602.00 E.9602.15 E.9600.01 E.9601.01 E.9602.01 E.9603.15

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Table of contents

Introduction	27
Legal Advices.....	27
Explanation of symbols.....	28
Security advices / warnings	29
Intended Use.....	32
Initial Operation.....	33
Technical Specification.....	34
Operation.....	36
Cleaning, disinfection, sterilization.....	37
Packaging.....	46
Trouble shooting.....	47
Warranty, service and repair.....	48
Disposal.....	49

Introduction

We thank you that you have decided to purchase nephroscopes of Gimmi® GmbH.

Endoscopes are medical devices that are manufactured to the highest technical standards and require careful handling, care and storage. Under the circumstances they will satisfy the high demands placed on them and will be usable to fulfil satisfaction over a long period time.

Gimmi® GmbH products are precision devices.

Please keep your endoscopes always with the utmost care, so you will enjoy them over a long period.

Read this manual carefully before using your new product. Thereby you keep yourself, the patient and any third parties from damage, which could occur by incorrect installation or by improper operation.

Legal Advices

This technical manual includes ownership protected information, which underlay copyright conditions. All rights are protected. This technical manual must not be copied, neither by photocopies, microfilms or other procedures and must not be distributed or saved, neither complete nor in extract and not without explicit written authorization of Gimmi® GmbH. Nomenclature which is at the same time a registered trademark is not especially marked. In the absence of trademark declaration, it may not be gathered that a designation will be a free trademark. Gimmi® GmbH would greatly appreciate being informed about any errors or omissions that may be found in the content of this instruction manual. Regarding continuous development and improvement of our products we reserve the right to alter technical features without written notice.

Explanation of symbols

Please take notice.

Warning Indicates to possible dangerous situations. No observance can lead to dead or serious injuries.

Attention Indicates to possible dangerous situations. No observance could cause injuries or product damages

Please consider the technical manual.
Please read technical manual before use.

Please consider the technical manual.
Please read technical manual before use.

Attention! / Caution!
Please consider the accompanying documents.

Reference number

Serial number

Date of manufacturing

Manufacturer

Non-Sterile – Sterilize before use

Temperature

Humidity

Pressure

Caution glass fragile

Do not expose to radiation

Do not expose to any weather

Devices have to be disposed separately from household/consumer waste

The recycling symbol serves to identify materials, which may be recycled due to the national regulation. The recycling code is built up by a circle of three arrows and the number of material inside the circle.

CE-Marking according to Council Directive 93/42/EEC.

CE-Marking

By CE-Marking the conformity of the medical device with the following directives is indicated:

Nephroscopes

are produced according to the requirements of:

DIN EN ISO 13485

Medical Devices

Council Directive 93/42/EEC according to Annex II

A circular CE marking with the number 0123 inside.**Security advices / warnings**

Improper use of this device and non-observance of the given instructions, warnings and precautions can lead to serious risks and consequences of the surgery or to injury, damage and even death of the patient, user and any third parties or damage to the endoscope!

This manual is designed to get familiar with the equipment and its intended capabilities in detail and use. The manual must therefore always be enclosed with the device.

The instruction manual contains important information needed for safe, proper, and efficient operation of the device.

We reserve the right to alter technical features, so that deviations of the content or pictorial representation are possible.

This manual should help to make the use of nephroscopes Gimmi® GmbH easier. But it is not designed as a guidance for endoscopic procedures and does not contain detailed description of endoscopy, and thus is not suitable for beginners to familiarize themselves with this surgical technique!

This device may only be used by technically competent and trained persons who are instructed to handle the instrument.

All persons using the devices have to read these instructions carefully first.

Handling of the device has to be carried out in accordance with these instructions.

Use the endoscope for endoscopic purposes only!

Inspect the endoscope, the appropriate endoscopic accessories and all devices connected to optical and mechanical parts, for all possible damages to exclude the risk of injury before every use! Defective and loose parts affect the function and safety and have to be removed immediately! Endoscopes with damaged or defective parts must not be used any longer. The full function and intended use of the medical device has to be ensured and verified by the user when using accessories and other components before use.

In case of doubt, contact your dealer or the manufacturer.

The device is not intended for use in potentially explosive areas!

Mechanical stress due to falling off, strong buckling, bending in a narrow winding, strong shocks, and torsion, tensile load or compressive stresses can result in damage or destruction of the endoscope and thus lead to malfunction!

We are not responsible for damages caused by misuse of the endoscope, no liability!

Endoscopic procedures should be performed by trained professionals (e.g. physicians) with appropriate training and experiences in performing endoscopic procedures only. It is within the responsibility of the user to consult information about indications, contraindications, potential complications, risks, and the development of endoscopic methods continuously.

Notes on combinations with other medical devices

The connection of other equipment or supplies (such as TV adapter, light sources, optical fibre cables, cameras, monitors, printers, video recorder, image processing systems, filing systems, pumps, shavers, insufflator, RF devices, work items, laser devices, pneumatic or electrohydraulic lithotripters etc.) opens up a variety of therapeutic applications. Follow the instructions and security advices of the used devices and accessories. Make sure that users are adequately trained. In case of doubt, contact your dealer or the manufacturer.

Protective measures for HF applications, including laser application, high-energy applications are not integrated in the Gimmi® GmbH device. Note that only the devices which are permitted for medical purposes may be adapted!

A thorough understanding of procedures used for endoscopic laser and electrosurgical treatment, applied principles and methods is needed to avoid shock and burn risks for patients and users as well as damage to other equipment and instruments. Liability claims arising from improper use or combination with other devices and instruments are excluded. Make sure when joint operation of an endoscope with electronic medical devices is performed, that the BF conditions (isolated, floating applied part) are observed.

If endoscopes should be used with electronic medical devices and / or energy-powered endoscopic usable accessories, leakage currents may be added up.

Notes for use with light sources

Gimmi® rigid endoscopes can be adapted to all common light sources for medical endoscopy. The malfunction of a used light source might lead to hazards. Keep an operational replacement light source available or use light sources that have a spare bulb. If bulb change is necessary during endoscopic application do not move the endoscope during the bulb change. Only if it is possible pull out the endoscope carefully for bulb change. Remember that light is an energy source that can heat each endoscope optics. The application time is limited by the selection of the light source.

In combination with high intensity light sources, both the light source side and the instrumental side optical fibre end can achieve temperatures that can cause burns. In addition, light of high energy radiation can lead to a temperature increase in tissue. During invasive application temperatures above +41 °C should be avoided, as this can cause tissue damage!

Therefore, avoid direct tissue contact and if applicable, pay attention to adequate irrigation of the operative field and the respective device-specific instructions and safety precautions.

Notes for use with high frequency surgical equipment

Prior to application of endoscopic high frequency treatment, surgical patients should be prepared in a suitable manner for the intended intervention. This includes activities to eliminate and to prevent the formation of ignitable gases in particular. In contrast to conventional high-frequency surgery inappropriate (particularly to low) power settings in high-frequency endoscopic surgery can cause a distinctive depth effect in the surrounding tissue.

The power adjustment should be made according to the users experience with respect to appropriate clinical references and / or appropriate training.

To avoid burns and / or unwanted depth effects in the surrounding tissue and to avoid endoscope damage, the high-frequency current should be switched on only if the appropriate application part (electrode) can be seen through the endoscope. The corresponding manuals, specifications and security advices should be respected. Never touch the endoscope while operating with an active electrode.

The insulation protection is provided by accessories and high frequency equipment.

**Notes for use with lasers**

If endoscopes or endoscopic accessories are used with laser devices, suitable protective glasses have to be worn to avoid potential damage to the eyes.

To avoid burns and / or unwanted depth effects in the surrounding tissue or damage of the endoscope, the laser power should be activated only if the tip of the laser fibre can be seen through the endoscope. The respective device-specific instructions and safety precautions have to be observed. Never touch the endoscope while operating with an activated laser fibre.

**Notes for use with lithotripters**

To avoid danger and in relation to possible restrictions on use of ultrasonic, electro-hydraulic, pneumatic, and mechanical lithotripsy device-specific instructions and safety precautions have to be observed. Appropriate surgery sheaths may be used for stone extraction with stone forceps. The required dimensions of suitable instrument can be gathered from technical specifications of the respective single devices. Never touch the endoscope while operating with an activated Lithotripsy unit.

**CT (Computer tomography)**

Certain metals of the endoscope can be dangerous due to heating during the application, so that an X-ray examination may be contraindicated in such patients. Due to X-radiation, optical components can discolour and thus can lead to endoscope damage. Concomitant use of CT (computed tomography) / X-ray and endoscopes can lead to hazards. Note therefore appropriate manufacturers and safety instructions.

**MRT (Magnetic resonance tomography)**

Due to magnetic field induced movements / relocations or heating some metals of the endoscope can be dangerous during the investigation, so that an MRI scan may be contraindicated in such patients. The optical and electrical medical devices for endoscopy may be damaged by magnets. Metals of the endoscope can cause side effects and visual disturbances. Concomitant use of MRI / magnetic resonance imaging and endoscopes can lead to hazards. Note therefore appropriate manufacturers and safety instructions.

Intended Use

1. General Description of the product

Rigid endoscopes of Gimmi® GmbH are optical devices which allow to inspect the insides of the body through natural body openings or artificial cavities and to visualize internal organs, tissues, and structures. To the human eye such insights remain hidden without the aid of an endoscope.

Endoscopes transmit the image information of the examined object or the examined orifice via a lens system inside the endoscope to the eyepiece. The light for the examination, inspection or operation is provided by a cold-light source, which can be connected to different standardized light source adapters. For the Gimmi® GmbH nephroscopes, we recommend light cables with a diameter of 4,8 mm. With the aid of fibre optics, it is possible to guide the light to the distal end of the endoscope without huge losses. The image at the eyepiece is captured from a camera and is displayed on a monitor.

Nephroscopes might be used for diagnostic purposes or in combination with surgical instruments for surgical purposes (minimally invasive). Therefore it's necessary to use a sheath (T.9600.30 or T.9600.33, depending on scope).

Normally, endoscopic procedures for diagnostic and therapeutic purpose are less stressful for the patient than traditional methods.

Endoscopes are intended for short-term use. Usually less than 24 hours.

For nephroscopes, the following parameters are variable:

- Diameter
- Working length
- Direction of view
- Field of view
- Diameter working channel
- Number of fluid channels

2. Medical Purpose

Endoscopes are intended for endoscopic visualization of internal anatomical structures of the patient. Nephroscopes are used to diagnose, treat or relieve renal disorders. The endoscopic shaft is inserted into the body percutaneously through a small skin incision (surgically invasive). The endoscope is gently pushed forward with the aid of an operation tube. This enables a view on the surgical area considering a maximum protection of the surrounding tissue structures.

3. Indications

Nephroscopes are used for example for the following minimally invasive surgical procedures:

Urology:

- Diagnostic percutaneous nephroscopy
- Therapeutic nephroscopy
- Percutaneous nephrolithotomy (PCNL)

4. Contraindications

The use of rigid Gimmi® GmbH endoscopes is contraindicated when endoscopic procedures are contraindicated.

5. Intended User

Nephroscopes must only be used by technically competent, medically qualified and in the functioning of the device – regarding minimally invasive surgery – trained medical staff. As well as staff of the processing and sterilization unit. Before using the nephroscope, the user has to read the corresponding instructions for use carefully.

6. Intended Patient

See point 3 and 4 for the corresponding indications and contraindications. Patients of all age needing endoscopic diagnosis and/or therapy at any area of the entire body where endoscopic diagnosis and therapy is required and manageable.

7. Intended Location of use

The usage environment corresponds to a clinic, e.g., an operating room or a processing and sterilization unit (CSSD). Pay attention to the interaction between technical devices. Network devices must not interfere or impair the functionality of other devices, especially in the field of medical technology.

Please read carefully the instruction manual prior to operation.


These instructions must be kept in a prominent position in proximity to the medical device.

Initial Operation

Remove the Gimmi® GmbH rigid endoscope and its accessories from the packaging and carefully remove all packing materials.

The delivery includes

1 piece	Rigid Endoscope
1 piece	Wolf-Adaptor
1 piece	Storz-Adaptor

The rigid endoscope and accessories must be checked for completeness and obvious damage immediately after unpacking. Damage can be claimed only if the supplier is notified immediately (within 24 hours).

Please use for any necessary return of equipment or accessories the original packaging. Describe the problem, identify the malfunction and designate a contact person for possible inquiries.

Caution


Please note the references made in section 5 (combinations with other medical devices) if equipment or other medical devices are connected to the endoscope.

Follow the instructions and safety advices of the used devices and the accessories strictly.


Gimmi® rigid endoscopes are supplied non-sterile and have to be cleaned before first use according to the processing instructions. Rigid endoscopes have to be disinfected or sterilized according to the medical indication (please note section "cleaning, disinfection, sterilization").

Technical Specification

1. Description of nephroscopes with working channel 45° version:




90° version:


Compatibility	45° Version	90° Version
Light port adaptor	90° angled 3 Adaptors according ISO/TS 18339:2015-11 Compatibility fixed ACMI; threaded Wolf and Storz/Olympus	90° angled 3 Adaptors according ISO/TS 18339:2015-11 Compatibility fixed ACMI; threaded Wolf and Storz/Olympus
Eyepiece	45° angled Ocular according ISO/TS 18339 threaded	90° double angled Ocular according ISO/TS 18339 threaded
Instrument port	Distal Shaft lock system, proximal M18x1, 0 thread	
Fluid port	2x Luer-Lock, according DIN EN ISO 80369-7	
Objective	CAD-designed tension-free composed-objective	
Image transmission	CAD-designed rod lens system for optimal sharpness, colour reproduction and resolution	
Light-guiding system	High quality light fibres	
Optical glass	Distal and proximal Sapphire	
Biocompatibility	All materials that come in contact with the body are biocompatible. All metal parts are made of high quality, medical grade stainless steel.	
Sterilization	Autoclavable, max. 138° C / 3 bar / 18 min	
Classification	Risk class IIa	
Certifications	Council Directive of 93/42/EEC, DIN EN ISO 13485	
CE-Marking	CE0123	

2. Definition of technical dimensions

<p>Diameter</p>	<p>Outer tube dimension</p> 	<p>in mm</p>
<p>Optic</p>		<p>in mm</p>
<p>Direction of view α</p>		<p>0° Colour code green</p> 
<p>Direction of view α'</p>		<p>15° Colour code white</p> 
<p>Field of view β</p>		<p>in degree</p>
<p>Working length</p>		<p>in mm</p>
<p>Complete length of working channel</p>		<p>in mm</p>
<p>Complete length</p>		<p>in mm</p>
<p>Complete height</p>		<p>in mm</p>

Operation


Setup of minimum configuration in endoscopy (this may vary due to different combination of units)

TV-Adapter connection to a camera

TV-Adapter connection to an endoscope


Pay attention to fibre optics adjustment!


Instrument connection

Pay attention to the correct positioning of the instrument to the optics!

Pay attention to instrument adjustment!

**Completely connected endoscope****Cleaning, disinfection, sterilization****1. General Principles / Introductory Remarks**

Also note that the manner of treatment can have significant impact on the life of endoscopes. Check your endoscope optics function (such as adequate lighting of the fibres, clear, sharp, bright and round picture) and check your endoscope for possible damage (such as sharp edges, loose parts or visible deformation of materials) before each use.

All Gimmi® GmbH endoscope optics have to be cleaned before each use, disinfected and sterilized, and this is especially true for first-time use after delivery, because all instruments are supplied non-sterile (cleaning and disinfection after removal of the protective packaging; sterilization after packaging). An effective cleaning and disinfection are essential prerequisites for effective sterilization.

As part of your responsibility for the sterility of the instruments / devices in use always ensure:

- that adequate equipment and product-specific validated procedures for cleaning / disinfection and sterilization are used in principle only,
- that used equipment (disinfector, sterilizer) is serviced and checked regularly and
- that validated parameters are strictly adhered to in each cycle.

During the application please pay special attention to the collection and separation of soiled / contaminated instruments. Keep them separated and do not place them back on the instrument tray in order to avoid higher contamination of the assembled instrument tray. Clean / disinfect contaminated instruments. Place them back on the instrument tray and sterilize the instrument tray fully stocked afterwards.

The endoscope optics should be cleaned immediately after each use. To avoid adverse effects on the components of the endoscope lens only use demineralised water for cleaning if possible. An effective cleaning and disinfection are essential prerequisites for effective sterilization. Before each use, the endoscope optics has to be cleaned, disinfected and sterilized. Please observe the applicable regulations of your country as well as the hygiene directives applicable to medical practices/hospitals.

This especially applies to the different guidelines / requirements for effective prion inactivation.

2. Instructions for processing

These instructions consider the procedure according to DIN EN ISO 17664 processing of medical devices (cleaning, disinfection and sterilization).

A washer-disinfector (WD) must be used for cleaning and disinfection. It is possible to use a manual cleaning as pre-cleaning process. But because of a significantly lower efficiency and reproducibility of the manual method, it is not allowed to apply a manual pre-cleaning without a subsequent mechanical cleaning and disinfection.

Before every cleaning, it is necessary to carry out a pre-treatment.

As endoscope optics are provided with irrigation and aspiration stopcocks, these should be disassembled for sterilization. We recommend a steam sterilization / autoclave process. When using a manual cleaning and disinfection procedure, product and process specific validation is the sole responsibility of the user.

3. Limitations and restrictions on processing, durability

Frequent processing of rigid endoscope optics has an impact on their usability. The end of product lifetime is usually determined by wear and damage caused by use.

Please note

With due care, rigid endoscopes can be used up to 100 times if they are undamaged and can still be cleaned. Any further re-use or the use of damaged and / or contaminated endoscopes is the responsibility of the user. However, if you are not sure whether your endoscope optics are operational, we recommend returning the rigid endoscope to our service address for review or replacement. After 100 sterilization cycles, the endoscopes must be sent to our service center for maintenance (see Chapter 13). Following the maintenance performed by us, the endoscope can be used again for 50 cycles.

4. Preparation at the site of use

In clinical practise used endoscope optics sometimes come in contact with corrosive etching agents and drugs. Endoscopes have to be removed from other equipment after use immediately. Surface contamination should be cleaned with a lint-free soft coarse cloth or paper towel.

Please note

"Dropping" of instruments (optics) during surgery will inevitably lead to damage.
Please pay special attention to place back the device gently after use.

5. Preparation before cleaning

In clinical practise used endoscope optics get sometimes in contact with corrosive etching agents and drugs. Coarse impurities on the instruments have to be removed directly after use (within 2 h). Immediately after the operation, a wet cleaning should be performed to prevent the blood, protein and other substances in the endoscope from drying out and to protect the personnel from these substances. Dried protein complicates the cleaning, disinfection and sterilization. Use running water or a disinfectant solution. The disinfectant should be aldehyde-free (otherwise there is the possibility of fixation of blood smears on the scope), have proven efficacy (eg. VAH / DGHM - or FDA approval or CE mark) and should be suitable for disinfecting instruments.

To remove impurities manually, use only a soft brush or a clean, soft cloth that is used only for this purpose. Never use metal brushes, steel wool or other sharp objects.

**Please note**

The disinfectant used in the pre-treatment is only for personnel's safety during this process. It does not replace the next disinfection step carried out after cleaning.

**Caution**




- Endoscopes must not be cleaned in an ultrasonic bath.
- Use approved cleaning agents only. Do not put the endoscope in alcohol or other corrosive liquids.

Procedure

Remove all sealing caps, valves, bushing grommets from the endoscope (discard disposable items immediately).

Remove all adaptors from the endoscope (discard disposable items immediately).

Rinse the endoscope with sufficient deionized water.

Please do not forget to re-mount the adapter and the flow valves after sterilization. Always grease the flow valves after sterilization. Only use greases which are approved for medical devices.

6. Cleaning and disinfection

Basics

A washer-disinfector (WD) must be used for cleaning and disinfection. It is possible to use a manual cleaning as pre-cleaning process. But because of a significantly lower efficiency and reproducibility of the manual method, it is not allowed to apply a manual pre-cleaning without a subsequent mechanical cleaning and disinfection.¹

Before every cleaning, it is necessary to carry out a pre-treatment.

When using a manual cleaning and disinfection procedure, product and process specific validation is the sole responsibility of the user.

Pre-Treatment

Coarse impurities on the instruments have to be removed directly after use (within 2 h). For this, directly after disconnection of the plug from the device, close the plug by the attached plug protective cap.

Procedure:

1. Check for closure of the plug with the plug protective cap.
2. Rinse the instruments at least 1 min under running water (temperature < 35 °C/95 °F). Push each button at least three times during pre-rinsing.
3. Soak the disassembled instruments for the given soaking time in the pre-cleaning solution², so that the instruments are sufficiently covered. Pay attention that there is no contact between the instruments or between instrument and plug. Assist cleaning by careful brushing of surfaces (but not the glass surface of the camera) with a soft brush (at beginning of soaking). Push each button at least three times during pre-cleaning.
4. Then, remove the instruments of the pre-cleaning solution and post-rinse them at least three times intensively (at least 1 min) with water. Push each button at least three times during post-rinsing

Pay attention to following points during selection of the cleaning detergent¹:

- basic suitability for cleaning instruments made of metal or plastic
- compatibility of the cleaning detergent with the instruments (see chapter „material resistance“)

Caution

Endoscopes must not be cleaned in an ultrasonic bath.

¹ When using a cleaning and disinfection detergent (e.g. in consequence of staff's safety) please consider that this should be aldehyde-free (otherwise there is the possibility of fixation of blood smears), should have proven efficacy (eg VAH / DGHM - or FDA approval or CE mark) and should be suitable for disinfecting instruments. Furthermore, it should be suitable for the disinfection of instruments made of metallic or plastic material and compatible with the instruments (see chapter "material resistance"). The disinfectant used in the pre-treatment is only for personnel's safety during this process. It does not replace the next disinfection step that has to be performed after cleaning.

Pay attention to the instructions of the detergent manufacturer regarding concentration, temperature and soaking time as well as post-rinsing. Only use freshly prepared solutions as well as sterile or low contaminated water (max. 10 germs/ml), low endotoxin contaminated water (max. 0.25 endotoxin units/ml), for example purified/highly purified water.

Additionally, use soft, clean, and lint-free cloths and/or filtered air for drying.

Manual cleaning and disinfection

When choosing the cleaning agent and disinfectant ensure:

- that these are basically suitable for the cleaning and disinfection of instruments made of metals and plastics,
- that a disinfectant with proven efficacy (e.g. VAH/DGHM- or FDA approval or CE mark) is used and that it is compatible with the cleaning agent,
- that chemicals used are compatible with the instruments (see "material resistance").

If possible combined cleaning / disinfecting agents should not be used. Combined cleaning / disinfecting agents can be used only in cases of very low contamination (no visible impurities).

Concentrations and contact times of cleaning and disinfectant agents specified by the manufacturer must be strictly adhered to.

Only use freshly prepared solutions and sterile or low contaminated water (max. 10 germs/ml), low endotoxin contaminated water (max. 0.25 endotoxin units/ml), e.g. purified water / highly purified water

Additionally, use filtered air for drying.

Procedure:

Cleaning:

1. Rinse endoscopes under running water. Remove surface contamination with a soft cloth.
2. Rinse the lumens five times with water using a disposable syringe (working channel 50 ml, rinsing and suction channel 10 ml).
3. Disassemble the endoscopes as far as possible. Remove the fiber optic connector, adapter, handle, irrigation valves, etc.
4. After pre-treatment, place the endoscopes in the cleaning solution for at least 5 minutes so that the instruments are sufficiently covered and thoroughly cleaned with a soft cloth or soft brush. Make sure that the instruments do not touch each other. Move the moving parts back and forth several times during cleaning.
5. Clean the lumens with a soft plastic brush.
6. If necessary, connect the rinsing adapter to the endoscopes.
7. At the beginning and at the end of the cleaning, rinse the lumens five times with cleaning solution (working channel 50 ml, rinsing and suction channel 10 ml).
8. Then rinse the products three times with sterile, deionized water.
9. Using a disposable syringe, rinse the lumens five times with sterile, deionized water (50 ml working channel, 10 ml irrigation and suction channel).
10. For cleaning the endoscopes, use a fresh cleaning solution.
11. Check the instruments (see chapter "Control" and "Maintenance").

Disinfection:

12. Disassemble the endoscopes as far as possible. Remove the fiber optic connector, adapter, handle, irrigation valves, etc.
13. If necessary, connect the rinsing adapter to the endoscopes.
14. Place the endoscopes in the disinfectant solution for 12 minutes so that the instruments are adequately covered. Make sure that the instruments do not touch each other.
15. At the beginning and at the end of the disinfection, rinse the lumens five times with disinfectant solution (working channel 50 ml, rinsing and suction channel 10 ml).
16. Then rinse endoscopes five times with sterile, deionized water.
17. Using a disposable syringe, rinse the lumens five times with sterile, deionized water (50 ml working channel, 10 ml irrigation and suction channel).
18. Use a fresh disinfectant solution for the disinfection of the endoscopes.
19. The evidence of the general suitability of the instruments/devices for effective manual cleaning and disinfection has been provided by an independent accredited test laboratory using the detergent Cidezime / Enzol and the disinfectant Cidex Opa (Johnson & Johnson GmbH, Norderstedt). The method described above has been considered.

Automated cleaning/disinfection (WD (Washer-Disinfector))

The automated reprocessing have been evaluated due to ISO 17664, AAMI TIR 12, AAMI TIR30, ASTM E 2314 mod. and ASTM E 1837 mod. These procedures are validated ² and suitable to reprocess the named endoscopes for another application. But be aware that it is still user's responsibility to ensure that actually applied reprocessing (equipment, materials and CSSD employees) achieves the desired outcome. To fulfil this criteria, a regular validation of the process is necessary. Gimmi® GmbH recommends to the user to comply with the named standards for reprocessing of medical endoscopes, that are listed above.

Pay attention to following points during selection of the WD:

- fundamentally approved efficiency of the WD (for example CE marking according to EN ISO 15883 or DGHM or FDA approval/clearance/registration)
- possibility for an approved program for thermal disinfection (A_0 value ≥ 3000 or – in case of older devices - at least 5 min at 90 °C/194 °F; in case of chemical disinfection danger of remnants of the disinfectant on the instruments)
- basic suitability of the program for instruments as well as sufficient rinsing steps in the program
- post-rinsing only with sterile or low contaminated water (max. 10 germs/ml, max. 0.25 endotoxin units/ml), for example purified/highly purified water
- use of filtered air (oil-free, low contamination with microorganisms and particles) for drying
- regularly maintenance and check/calibration of the WD

Pay attention to following points during selection of the cleaning detergent:

- fundamental suitability for the cleaning of instruments made of metallic or plastic material
- additional use - if thermal disinfection is not used - of a suitable disinfectant with approved efficacy (for example VAH/DGHM or FDA/EPA approval/clearance/registration or CE mark) compatible to the used cleaning detergent
- compatibility of the used detergents with the instruments (see chapter „material resistance.“)

Pay attention to the instructions of the detergent manufacturers regarding concentration, temperature and soaking time as well as post-rinsing.

Procedure:

1. Disassemble the endoscopes as far as possible. Remove the fiber optic connector, adapter, handle, irrigation valves, etc.
2. If necessary, connect the rinsing adapter to the endoscopes.
3. Then position the pretreated endoscopes in a washer-disinfector (Disinfector G 7836 CD, Miele, Gütersloh). Make sure that the instruments do not touch each other.
4. Connect the lumens with suitable, flexible irrigation tubing and adapters.
5. Start the program.
6. The machine cleaning was carried out using "neodisher MediClean forte" as a cleaning agent.
7. The machine is cleaned at 50 ± 2 °C for at least 5 min. The thermal disinfection was carried out at 90 ± 2 °C for at least 5 min.
8. After completion of the machine cleaning / disinfection, remove the endoscopes from the washer-disinfector under low-germ conditions.
9. Note:
The fundamental suitability of the instruments for an effective automated cleaning and disinfection was demonstrated by an independent, governmentally accredited and recognized (§ 15 (5) MPG) test laboratory by application of the WD G 7836 CD, Miele & Cie. GmbH & Co., Gütersloh, (thermal disinfection) and the pre-cleaning and cleaning detergent Neodisher medizym (Dr. Weigert GmbH & Co. KG, Hamburg) considering to the specified procedure.

² The reprocessing validation was performed by Medical Device Services - Dr. Rossberger GmbH in 82205 Gilching. It can be viewed at any time under document numbers 191957-10 and 191958-10.

Drying

The endoscope and the camera head must be completely dry after disinfecting. Disinfected endoscopes should always be stored in a closed sterile area or cabinet and protected against heat, radiation, dust, moisture, temperature changes and contamination.

Check

Check the device after cleaning or cleaning / disinfection for corrosion, damaged surfaces, impurities and functionality, respectively. Do not use a damaged device (limiting of the number of re-use cycles see chapter "reusability"). If the device is still dirty, it must be cleaned and disinfected again.

Inspection of the mechanics and endoscope surface

The endoscope surfaces have to be undamaged and in particular free of sharp edges. Check for dents, bends, mechanical / thermal damage caused by radiofrequency or laser surgery equipment as well as for cracks and spalling.

Inspection of the fibre optics

Hold the distal endoscope end toward a lighted window or a bright ceiling light.

Caution

Do not use a cold-light source for this test. Direct view into the radiated light from a cold-light source can cause eye damage.

Look at the light guide connector. The individual fibres now appear bright. Move the side, which is facing the bright ceiling light, slightly up and down. The brightness of the fibres changes a bit. It is uncritical if individual fibres remain dark. A breakage rate of about 20 to 30% impedes the endoscopic procedures severely.

The surfaces of the light entry and exit surfaces should be smooth and clean. Rough surfaces with deposits, tangible or withdrawn individual fibres may lead to insufficient lighting. Further use and processing may result in progressive damage to the endoscope.

Caution

Endoscopes with damaged fibre optics should be returned to the manufacturer or an authorized service specialist for checking.

Verification of the proximal and distal areas of glass

Glass surfaces have to be clean and debris-free. Persistent encrustation, observed during visual examination, should be removed with appropriate cleaning pastes or alcohol-soaked cotton swab or toothpick. Inadequate rinsing of the optics after cleaning and disinfection is often the cause of precipitates.

Corresponding to the indications appropriate working distance the image has to be sharp and clear. A fuzzy, non-circular, cloudy or foggy image points out to a damage of the scope.

Caution

- Endoscopes with indelible persistent encrustations should be returned to the manufacturer or an authorized service specialist for checking.
- Endoscopes with damaged glass surfaces (e.g. chips), impaired image quality or noticeable surface damage and deformations may no longer be used. They should be discarded or sent back to the manufacturer or an authorized service specialist for checking.

Maintenance

Instrument oils or grease must not be used. Put the disassembled endoscopes (fibre optic light guide, handles) together again.

Please note

Endoscopes do not require regular maintenance carried out by the manufacturer.

Packaging

Please insert the cleaned and disinfected device in single-use sterilization packaging's (single or double packaging), which fulfil the following requirements (material/process):

- ANSI AAMI ISO 11607 (for USA: FDA clearance)
- suitable for steam sterilization (temperature resistance up to at least 138 °C (280 °F), sufficient steam permeability)
- sufficient protection of the instruments as well as of the sterilization packaging against mechanical damage

Pay attention that there is no contact between the device or between device and plug.
Avoid contact of the aluminium ferrule with the plug and camera surface.
Avoid scratches on the surface.

Sterilization

Please use only the listed sterilization procedures for sterilization; other sterilization procedures must not be used.

Steam sterilization:

- fractionated vacuum/dynamic air removal procedure^{3 4}(with sufficient product drying⁵)
- steam sterilizer according to EN 13060/EN 285 or ANSI AAMI ST79 (for USA: FDA clearance)
- validated according to EN ISO 17665 (valid IQ/OQ (commissioning) and product specific performance qualification (PQ))
- maximum sterilization temperature of 134 °C (273 °F) plus tolerance according to EN ISO 17665
- sterilization time (exposure time at the sterilization temperature):

Area	fractionated vacuum/dynamic air removal	gravity displacement
USA	at least 4 min at 132 °C (270 °F), drying time at least 20 min	not recommended
other countries	at least 5 min ⁶ at 132 °C (270 °F) / 134 °C (273 °F)	not recommended

The flash/immediate use sterilization procedure must not be used.

Do not use dry heat sterilization, radiation sterilization, formaldehyde and ethylene oxide sterilization, as well as plasma sterilization.

This results in the following parameters for use on our endoscopes:
max.138°C

AUTOCLAVE

max.3 bar / 18 min

The use of less effective gravitational method is permitted for non-availability of the fractionated vacuum method only, can have significantly longer exposure times and has to be validated by the user as to suitability and effectiveness.

The evidence of the general suitability of the instruments/devices for effective steam sterilization has been provided by an independent accredited test laboratory using the steam sterilizer Systec V-150 (System GmbH Labor Systemtechnik, Wettenberg) and the fractionated vacuum procedure. At this, typical conditions in clinics and medical practices as well as the method described above has been considered.

3 at least three vacuum steps

4 The less effective gravity displacement procedure must not be used in case of availability of the fractionated vacuum procedure, requires significantly longer sterilization times as well as a sterilizer, procedure, parameter, and product specific validation under sole responsibility of the user.

5 The effectively required drying time depends directly on parameters in sole responsibility of the user (load configuration and density, sterilizer conditions, ...) and is thus to be determined by the user. Nevertheless, drying times less than 20 min must not be applied.

6 respectively 18 min (inactivation of prions, not relevant for USA)

**Caution**

- Flash sterilization is not permitted.
- Do not apply hot air sterilization, radiation sterilization, formaldehyde or ethylene oxide, or plasma sterilization.
- Other settings of autoclave and autoclave cycles might have negative effects on the device or its components.
- Weight and load of the sterilization material must not be exceeded as it may cause excessive condensate and result in rust damage.
- When sterilizing multiple instruments in a sterilization cycle, do not exceed the maximum load of the sterilizer.
- Please note that with increasing number of instruments in a sterilization cycle the success of sterilization will be decreased. Please refer to the manufacturer's instructions.
- The dryness of the optics must be achieved after cooling to room temperature.
- Please follow the sterilizer manufacturer's instructions, especially the ventilation times after sterilization. The relevant national legal regulations have to be observed.

Hint:

It is the user's responsibility to implement the listed sterilization processes in order to achieve the desired and required sterilization effects.

The instructions of the sterilizer manufacturer have to be strictly adhered to

**Note**

Please do not forget to re-attach the adapter, the one-way cocks and other accessories after sterilization. The one-way cocks have to be greased always after sterilization. Only use approved medical sterile greases.

Storage

Please store the device after sterilization in the sterilization packaging at a dry and dust-free place.

Material resistance

Be careful by selecting cleaning and disinfecting agents. Ensure that the following substances are not included:

- organic, mineral and oxidizing acids (minimum admissible value pH 5,5)
- strong alkalis (maximum admissible value pH 11, neutral / enzymatic or slightly alkaline detergent recommended)
- organic solvents (for example: alcohols, ethers, ketones, benzenes)
- oxidizing agents (for example: hydrogen peroxide)
- halogens (chlorine, iodine, bromine)
- aromatic / halogenated hydrocarbons
- oils

When selecting detergents, please keep in mind that corrosion inhibitors, neutralizing agents and / or rinse aids may cause critical remnants on the instruments. Acid neutralizing agents or rinse aids must not be applied.

Please do not clean any instruments by use of metal brushes or steel wool.

Please do not expose any instruments to temperatures higher than 138 °C (280 °F).

Reusability

Frequent processing of rigid endoscope optics has an impact on their usability. The end of product lifetime is usually determined by wear and damage caused by use.

**Please note**

With due care, rigid endoscopes can be used up to 100 times if they are undamaged and can still be cleaned. Any further re-use or the use of damaged and / or contaminated endoscopes is the responsibility of the user. However, if you are not sure whether your endoscope optics are operational, we recommend returning the rigid endoscope to our service address for review or replacement. After 100 sterilization cycles, the endoscopes must be sent to our service center for maintenance (see Chapter 13). Following the maintenance performed by us, the endoscope can be used again for 50 cycles.

Packaging


Caution

Gimmi® GmbH rigid endoscopes are supplied non-sterile

The original package includes

1 piece	Cover topside made of cardboard
1 piece	Damping-top part made of foam
1 piece	Damping-bottom made of foam
1 piece	Bottom-lower shell made of cardboard

Please use the original packaging for the possible return of endoscopes or accessories. Make sure that only sterile endoscopes are returned.

Always store non-sterilized endoscopes protected from heat, radiation, dust, moisture and temperature fluctuations. Always store sterilized endoscopes in sterile and suitable containers.

Store closed containers in a sterile area or cabinet protected from heat, radiation, dust, moisture, temperature fluctuations and contamination.

Storage

- Store non-sterilized endoscopes always protected from heat, radiation, dust, moisture and temperature changes.
- Sterilized endoscopes should always be stored in sterile containers or storage containers in a closed sterile area or cabinet safe and protected from heat, radiation, dust, moisture and temperature fluctuations.
- Temperature range of the storage should be between 5°C till 40°C.
- The humidity should be between 10% till 90%.
- The air pressure should be between 50kPa till 107kPa.

Transport

- When returning endoscopes or accessories, please always use the original packaging.
- Please mark the packaging with: !Caution Risk of breakage!
- Avoid strong vibrations and shocks.
- Avoid piercing objects through the box.
- Make sure that only sterile endoscopes are returned.
- Temperature range of the storage should be between -10°C till 60°C.
- The humidity should be between 5% till 95%.
- The air pressure should be between 50kPa till 107kPa.

Trouble shooting

Problem	Possible cause	Remedying of defect
Picture cloudy, foggy	Glass surfaces contaminated	Cleaning of glass surfaces according to section "manual cleaning"
	Deposits, coarse encrustations of glass surfaces	Remove deposits according to section "check water quality"
	Leaky, defective lens system	Send in the endoscope for repair
Picture too dark, too small illumination	Glass surfaces contaminated	Cleaning of glass surfaces according to section "manual cleaning"
	Deposits, coarse encrustations of glass surfaces	Remove deposits according to section "check water quality"
	Wrong light conducting cable connector	Check light conducting cable connector, replace if necessary
	Fibre optics defect	Check fibre optics according to section "Check"
	Defect light conducting cable, light source	Check light conducting cable, light source
Yellowish lighting	Dirty fibre optics	Cleaning of glass surfaces according to section "manual cleaning". If necessary, send in the endoscope for service.
	Dirty, broken light conducting cable	Check light conducting cable (for example, shine on white surface), replace if necessary
Staining, discoloration	Inadequate cleaning (for example, remaining protein residues)	Clean up, possibly with thorough scrubbing
	Inadequate rinsing of endoscope between treatment phases (especially before sterilization)	Ensure thorough rinsing between the treatment phases (see section "cleaning and disinfection" and "Drying"
	Contaminated, too often used disinfectants and cleaning solutions	Replace disinfection and cleaning solutions regularly
Leakage	Leaking connections	Check connections between sealing cap and irrigation stopcock
	Defect irrigation stopcocks	Send in the endoscope for repair

Warranty, service and repair

The Gimmi® GmbH provides 12 months warranty of the rigid endoscope. The duration of this warranty is limited to claims that are submitted within the specified warranty period from date of purchase of the endoscope, possibly related to repairs, stating the invoice number. This warranty applies to defects only that are not normal wear and tear, misuse, mishandling, improper or inadequate treatment or due to force majeure. In cases of maintenance or repair, please contact the Gimmi® GmbH service or an authorized repair specialist:

Gimmi® GmbH

Carl-Zeiss-Str. 6
78532 Tuttlingen / Germany
Tel.: +49 (0)7461 965 90 - 0
Fax: +49 (0)7461 965 90 - 33
Email: info@gimmi.de
www.gimmi.de

In the interest of expediting your request, please send the product with the following information:

- Part number (REF)
- Serial number (SN)
- Detailed error report

Specialty Retailer:

Manufacturer:

Gimmi® GmbH
Carl-Zeiss-Str. 6
78532 Tuttlingen / Germany
Tel.: +49 (0)7461 965 90-0
Fax.: +49 (0)7461 965 90-33
E-mail: info@gimmi.de
www.gimmi.de

If you have any further questions, please feel free to contact us. We are happy to help!

Disposal

Be aware of sharp edges during disposal to prevent injury.

The transposition of European legislation into national law, requires proper disposal of medical devices and instruments. The symbol below indicates that medical products must be disposed of separately from household and residual waste. Please dispose of medical products according to local and national laws.

Dispose of contaminated medical equipment at an appropriate collection point.

Medical waste is classified according to the Dangerous Goods Regulation for UN number "UN 3291" (Medical Waste)).

Please mark the waste container with the symbol below.

Dispose of medical waste according to local and national law.

You can dispose of uncontaminated packaging material according to local and national recycling regulations.

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