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Content

Symbols	15
Product Description / Intended Use	16
Intended User	16
Contraindication	16
Warning.....	17
Caution	17
Assembling The Gimmi® AlphaQuick System with Peek Handle Without Lock	18
Disassembling The Gimmi® AlphaQuick System with Peek Handle Without Lock	19
Assembling The Gimmi® AlphaQuick System with Peek Multifunctional Handle	19
Working Principle of The Peek Multifunctional Handle	20
Disassembling The Gimmi® AlphaQuick System with Peek Multifunctional Handle	21
Maintenance Instructions / Disinfection / Sterilization	21
Reprocessing Instructions for Automated Cleaning Process	22
Additional Instructions	24
Disposal.....	24
Guarantee and Warranty.....	24

Laparoscopic Instruments

Handling Instructions

Please read before use!

Symbols

-  Please do read the manual.

-  Precautions and Warnings

-  The instrument will be delivered unsterile and has to be prepared before each use according to our manual.

-  Storage conditions

-  Warning!
Instrument is sharp/pointed at the working end

Attention!
Instrument can be used with HF-electricity

Attention!
Wrong handling can cause a flashover of the HF-electricity.

Wear suitable personal protective equipment: gloves, goggles, etc.

Product can be disposed in the hospital's inhouse disposal system.

Federal law restricts this device to sale by or on order of a physician!

Product Description / Intended Use

Instruments for laparoscopic surgery are intended for cutting, preparing, resecting, gripping and sewing of body tissue. Read the following instructions in the sequence given in order to achieve optimum results in maintenance, cleaning and sterilization. Prior to using the instruments observe all Warnings, Cautions and Notes.

Patient population

Patients	Human medicine
Age	newborn to geriatric patients
Weight	> 2,5 kg
Health condition	Relevant, because patient is treated under anesthesia.

Intended User

The products must be used only in medical facilities by trained and skilled medical personnel. The products must not be used if, according to a qualified physician, the general condition of the patient is not adequate or if the endoscopic methods are contraindicated.

Contraindication

The using of instrument is contraindicated:

- when the therapeutic procedure is not accepted.
- when the intended use is differing from the real use.
- for patients with pacemaker.

Do not use the instrument if, in the opinion of the attending physician, the risks to the patient outweigh the benefits. The instrument is also contraindicated for cardiovascular system and central nervous system.

Warning

Immediately after unpacking inspect the instruments for transport damage and check their proper function. Any damage must be reported without delay.

The product is shipped in non-sterile condition. Therefore, instruments must be cleaned, disinfected and sterilized before each use. The jaw parts may be sharp or pointed – be careful when handling. When using laparoscopic instruments together with an HF generator always test the HF output power on the basis of the user's experience, clinical references and/or relevant training and education.

When using HF equipment always follow the HF unit manufacturer's instructions!
Maximum output voltage of the generator, U_{max} : 1300 Vp

The instruments must not be used under the following conditions:

- if the electrode/instrument insert is loose or damaged, the insulation of the shaft shows cracks or the connector is loose or wrongly connected. The use of instruments with damaged insulation, defective electrodes/instrument inserts or defective connectors may cause uncontrolled electro-surgical burns.
- if scratches, surface defects or fissures are visible on the instrument housing. The use of a defective instrument or an instrument that becomes defective during surgery may result in loose parts that could drop into the field of operation.
- if the connecting cable does not completely cover the electrical connector. This may result in an uncontrolled current flow.
- if large quantities of residue build up at the tool portion or if explosive gas concentrations are present in the field of operation.
- Current-carrying parts must not be touched

Always check the compatibility of the instruments and ensure an uninterrupted equipment ground connection before using the system.

Caution

In order to prevent the instrument from damage it should be disassembled and sterilized in open position.

Quick sterilization methods **must not** be used as the standard sterilization procedure since this may cause damage to the instrument.

Do not use excessive force when handling the instrument since this may cause the jaw parts or the instrument to break.

Applying undo force to turn the instrument while the tool portion is under load may damage the instrument insert.

Excessive force on the handles while the instrument insert is holding on to an object may damage the insert.

The instrument functionality should be checked before the surgeon uses the instrument. Should there be any indication of a malfunction the instrument must not be used and should be returned immediately to Gimmi® for checking and repair.

Repairs of Gimmi® products by non-authorized personnel can lead to malfunctions of the instrument and/or breakage of the jaw parts.

Expected product life

Designation	Value
Conditioning	≤ 50 cycles
Time	≤ 2 years


Assembling The Gimmi® AlphaQuick System with Peek Handle Without Lock


Press down the safety pin (6) on the moving arm (5) and **simultaneously fully open the handle arms.**

This is the correct position for attaching an instrument insert.

With the jaw parts closed slide the instrument insert into the shaft all the way to the stop **while keeping the handle arms fully open.** The round pins on the instrument insert (14) must engage in the recesses (13) of the shaft. Only hold the non-moving arm (4) in position. The moving arm (5) closes when the instrument insert hits the stop.

Close the instrument handle by applying some force to arm (5) until the instrument insert engages audibly.


After correctly positioning and engaging the insert the instrument is ready for use. Confirm the proper instrument function.


When using HF equipment always follow the HF unit manufacturer's instructions regarding output power and equipment connections.

Disassembling The Gimmi® AlphaQuick System with Peek Handle Without Lock


Press down the safety pin (6) at the moving handle arm (5) and **simultaneously fully open the instrument handle** (by applying some force) in order to release the instrument insert from its engaging mechanism.

Pull out the instrument insert (12).

Assembling The Gimmi® AlphaQuick System with Peek Multifunctional Handle


Deactivate the lock (3) by pressing the lever (1) in the direction indicated by the arrow until it engages. Press down the safety pin (6) on the moving arm (5) and **simultaneously fully open the instrument handle** (by applying some force).

This is the correct position for inserting an instrument insert.

With the jaw parts closed slide the instrument insert into the shaft all the way to the stop **while keeping the instrument handle fully open**. The round pins on the tang (14) must engage in the recesses (13) of the shaft. Only hold the non-moving handle arm (4) in position. The moving arm (5) closes when the instrument insert hits the stop.


Close the instrument handle by applying some force to handle arm (5) until the instrument insert engages audibly. After correctly positioning and engaging the insert the instrument is ready for use. Confirm the proper instrument function.


When using HF equipment always follow the HF unit manufacturer's instructions regarding output power and equipment connections

Working Principle of The Peek Multifunctional Handle

For activating the lock (3) press down the lock actuator (2) and the lever (1) will move automatically outward.

In order to release the lock press down the lock actuator (2).


For deactivating the lock (3) permanently simply push the lever (1) in the direction indicated by the arrow until it engages.


Note:

For sterilization always activate the lock

Disassembling The Gimmi® AlphaQuick System with Peek Multifunctional Handle

Deactivate the lock (3) by pressing the lever (1) in the direction indicated by the arrow until it engages. Press down the safety pin (6) at the moving handle arm (5) and **simultaneously fully open the instrument handle** (by applying some force) in order to release the instrument insert from its engaging mechanism.

Pull out the instrument insert (12).

Maintenance Instructions / Disinfection / Sterilization

Warning

1. Wear suitable protective clothing (gloves, eye protection, etc.) when using, cleaning or sterilizing this product.
2. Hand-held laparoscopic instruments from Gimmi® are shipped in **non-sterile**-condition. Therefore, before the first and any subsequent use clean, disinfect and sterilize the instrument according to the instructions below.
3. During use the instrument must not contact other metallic parts. Any contact with metal parts may damage the instrument. The use of defective instruments can result in loose parts that could drop into the field of operation.
4. Before returning an instrument for maintenance or warranty work it must be cleaned, disinfected and sterilized.

Caution

1. The quality of the water used for cleaning and steam sterilization influences the performance and the service life of the instruments. Only use purified water.

2. Do not use brushes or cloth containing metallic or abrasive substances for cleaning and disinfection of instruments since this may cause permanent scratches or damages.
3. Always observe the time of exposure to solvents specified by the solvent manufacturer in order to prevent corrosion and other damages.
4. The instruments must not be immersed in saline as this may cause corrosion.
5. Do not use strong acidic or basic solutions as these can cause damage to the instrument. For use of alkaline solutions, do not exceed pH=12.
6. When cleaning instruments in an ultrasound unit do not put them on the chamber bottom since this may cause damage or incomplete cleaning. Always follow the instructions of the equipment manufacturer.
7. When disinfecting or sterilizing instruments with multifunctional handle always activate the lock.

Note

1. The service life of the product depends primarily on its wear, the damages due to rough handling and the adherence to the intended use.
2. Damages caused by incorrect handling are not covered by warranty.

Reprocessing Instructions for Automated Cleaning Process

Advice

Reprocessing procedures have only limited implications to a surgical instrument. The limitation of the numbers of reprocessing procedures is therefore determined by the function / wear of the device. Used devices need to be reprocessed before sending back to the manufacturer.

Preparation at the point of use

Remove excess soiling by submerging the instrument in cold water (<40°C) immediately after use. Do not use a fixating detergent or hot water (>40°C) as this can cause the fixation of residues which may influence the result of the reprocessing process.

Transportation

Safe storage and transportation, for example in a closed container, to the reprocessing area to avoid any damage on the products and contamination to the environment.

Preparation for decontamination/ Manual pre-cleaning

Wear suitable protective clothing: gloves, eye protection, etc. The instruments must be reprocessed in an opened or if possible in a disassembled state. Immerse the instrument in cold tap water for at least 5 minutes. Brush (do not use a metal brush) the instrument under cold tap water until all visible residues are removed. Inner lumens, threads and holes are individually flushed with a water jet pistol for minimum 10 seconds in the pulsed mode. We recommend our T.1003.00 cleaning system.

Ultrasonic cleaning

Immerse the product in an ultrasonic bath with alkaline or enzymatic detergent (0,5%) and treat with ultrasonic for 15 minutes at 40°C (35 – 45 KHz). Turn the instrument several times during this time. Remove the product out of the bath and rinse with cold tap water. Follow the information provided by the equipment manufacturer.

Automated Cleaning

To optimize the cleaning process we recommend using of distilled, demineralized or purified water. Open instruments and carefully place them into a suitable screen basket.

Take-apart instruments must be taken apart for cleaning. If possible, connect instruments or instrument parts with lumens or channels using an existing Luer-Lock connection to the wash-up system of the cleaning and disinfection unit.

1. Rinse 1 minute with cold water
2. Drain
3. Rinse 3 minutes with cold water

4. Drain
5. Clean 5 minutes with 0.5 % alkaline detergent at 55°C or with enzymatic detergent at 45°C
6. Drain
7. Neutralize for 3 min with warm water (>40°C) and neutralizer
8. Drain
9. Rinse 2 minutes with warm water (>40°C)
10. Drain

Follow the operation instructions of the equipment and detergent's manufacturer.

Disinfection

Run automated thermal disinfection cycle in washer/disinfector as regulated by the national requirements and specified by Ao-Value (see ISO 15883).

Rinse with distilled, demineralized or purified water. Follow the information provided by the equipment manufacturer.

Drying

Drying of instrument through drying cycle of washer/disinfector. If needed, additional manual drying can be performed through lint free towel. Insufflate cavities of products by using residue-free compressed air. Follow the information provided by the equipment manufacturer.

Functional Test, Maintenance and Care

Let instruments cool down to room temperature and make sure that all joints and sliding surfaces are dry. Inspect instruments visually after each cleaning and disinfection for cleanliness. If necessary, repeat cleaning and disinfection. In case of electro-surgical instruments check the electrical insulation of the insert and/or the handle for burns and mechanical damage. Only use completely dry instruments. Instruments, such as clamps, scissors, needle holders, etc. must be specifically treated at the joints, locks, cutting edges and sliding surfaces. Oiling the entire product is unnecessary. Excess oil should be removed.

Oiling provides the product with optimal lubrication and protects against corrosion.

Gimmi® Instrument Oil J.8860.08 is suitable for steam sterilization and dry heat sterilization for up to 180°C.

Check for compatibility with corresponding parts. Inspect products for damages, loose parts, and wear. Perform a function test according to the instructions of use.

Sort out damaged instruments.

Packaging

Appropriate packaging for sterilization according ISO 11607 and EN 868.

Sterilization

Sterilization of products by applying a pre-vacuum process (according ISO 17665) under consideration of the respective country requirements.

Parameters for the pre-vacuum cycle:

- 3 prevacuum phases with at least 60 millibar.
- Heat up to a minimum sterilization temperature of 134°C; maximum temperature 137°C
- Pressure: 2.0 - 2.3 bar
- Minimum holding time: 4 min
- Drying time: minimum 10 min

The AlphaQuick System can be autoclaved in assembled condition by open flushing port – if the described cleaning process was used.

Follow the information provided by the equipment manufacturer.

Storage

Storage of sterilized instruments in appropriate containers, trays, etc. in a dry, clean and dust free environment at temperatures of 5°C to 40°C.

Reprocessing validation study information

The following test devices, materials and equipment have been used in this validation study;

- **Detergent:**
 - Neodisher FA; Dr. Weigert; Hamburg (alkaline)
 - Endozime, Fa. Ruhof (Enzymatic)
- **Neutralizer:**
 - Neodisher Z; Dr. Weigert, Hamburg
- **Washer/ Disinfector:**
 - Miele 7735 CD
 - Instrument Rack Miele E 327-06
 - Key Hole Surgery Rack Miele E 450

Details see validation reports

- **Cleaning:**
 - 01707011901-2/ 01707011901-3
- **Sterilization:**
 - 17308022611/ 09m111

Additional Instructions

If the described chemicals and equipment are not available, it is the duty of the user to validate his process.

The user must ensure that the reprocessing process, including resources, material and personnel, is suitable to achieve the necessary results.

The state of technology and national laws require that validated processes are complied with.

Disposal

This product can be disposed off in the hospital's in-house disposal system and can therefore be reintroduced to the raw material cycle. All national regulations on disposal must be observed.

Guarantee and Warranty

Gimmi® GmbH exclusively supplies tested and faultless products to its customers. All products are designed and manufactured to comply with maximum quality requirements. We refuse any liability for products which have been modified as compared to the original product, misused or handled or used improperly.

In case of repairs by companies, which are not authorised for repairs by Gimmi®, the warranty is inapplicable.

If a serious incident occurs in connection with a product of the company Gimmi®, the manufacturer must be notified immediately. After consultation with the latter, a notification is sent to the competent authority of the member state in which the user is established.

Gimmi® GmbH guarantees that this product is free of defects in material and design at the time of purchase. The warranty period is 2 years from date of purchase.

This guarantee applies to all purchases and is limited to free-of-charge repair or replacement of the product if it has been established that a material or manufacturing fault existed.

Normal wear and tear damage or misuse in the intended purposes are exempt from guarantee.

Shipping costs and risk of loss or damages are not covered by Gimmi® GmbH.

Prior to returning the product for maintenance purposes or for warranty claims, the product must be cleaned, disinfected, and sterilized.

Disclaimer

The manufacturer does not assume liability for immediate or secondary damages and the warranty claim is void if

- the product was improperly used, reprocessed or serviced,
- instructions and provisions stated in the operating manual were not observed,
- non-authorized personnel performed repairs, calibrations or modified the product in any manner.

The company Gimmi® assumes no liability, if this customer information is demonstrably violated.

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