



EN Instructions for use HF forceps



REF

Bipolar Flat Plug

B.9100.00	B.9100.18	B.9100.36	B.9100.54	B.9100.72	B.9100.90	B.9101.06	B.9101.24	B.9101.40
B.9100.02	B.9100.20	B.9100.38	B.9100.56	B.9100.74	B.9100.92	B.9101.08	B.9101.26	B.9101.42
B.9100.04	B.9100.22	B.9100.40	B.9100.58	B.9100.76	B.9100.94	B.9101.10	B.9101.28	B.9101.44
B.9100.06	B.9100.24	B.9100.42	B.9100.60	B.9100.78	B.9100.96	B.9101.12	B.9101.30	B.9101.46
B.9100.08	B.9100.26	B.9100.44	B.9100.62	B.9100.80	B.9100.98	B.9101.14	B.9101.32	B.9101.62
B.9100.10	B.9100.28	B.9100.46	B.9100.64	B.9100.82	B.9100.99	B.9101.16	B.9101.34	B.9101.64
B.9100.12	B.9100.30	B.9100.48	B.9100.66	B.9100.84	B.9101.00	B.9101.18	B.9101.36	B.9101.66
B.9100.14	B.9100.32	B.9100.50	B.9100.68	B.9100.86	B.9101.02	B.9101.20	B.9101.38	B.9101.68
B.9100.16	B.9100.34	B.9100.52	B.9100.70	B.9100.88	B.9101.04	B.9101.22		

Bipolar US Pin

B.9199.00	B.9199.18	B.9199.36	B.9199.54	B.9199.72	B.9199.90	B.9200.06	B.9200.24	B.9201.40
B.9199.02	B.9199.20	B.9199.38	B.9199.56	B.9199.74	B.9199.92	B.9200.08	B.9200.26	B.9201.42
B.9199.04	B.9199.22	B.9199.40	B.9199.58	B.9199.76	B.9199.94	B.9200.10	B.9200.28	B.9201.44
B.9199.06	B.9199.24	B.9199.42	B.9199.60	B.9199.78	B.9199.96	B.9200.12	B.9200.30	B.9201.46
B.9199.08	B.9199.26	B.9199.44	B.9199.62	B.9199.80	B.9199.98	B.9200.14	B.9200.32	B.9201.62
B.9199.10	B.9199.28	B.9199.46	B.9199.64	B.9199.82	B.9199.99	B.9200.16	B.9200.34	B.9201.64
B.9199.12	B.9199.30	B.9199.48	B.9199.66	B.9199.84	B.9200.00	B.9200.18	B.9201.36	B.9201.66
B.9199.14	B.9199.32	B.9199.50	B.9199.68	B.9199.86	B.9200.02	B.9200.20	B.9201.38	B.9201.68
B.9199.16	B.9199.34	B.9199.52	B.9199.70	B.9199.88	B.9200.04	B.9200.22		

Bipolar Non-Stick Flat Plug

B.9120.00	B.9120.08	B.9120.18	B.9120.32	B.9120.66	B.9120.74	B.9120.82	B.9120.98	B.9121.38
B.9120.01	B.9120.10	B.9120.19	B.9120.34	B.9120.67	B.9120.76	B.9120.84	B.9121.02	B.9121.39
B.9120.02	B.9120.12	B.9120.20	B.9120.36	B.9120.68	B.9120.78	B.9120.85	B.9121.20	B.9121.40
B.9120.04	B.9120.13	B.9120.22	B.9120.37	B.9120.70	B.9120.79	B.9120.86	B.9121.26	B.9121.42
B.9120.06	B.9120.14	B.9120.30	B.9120.38	B.9120.72	B.9120.80	B.9120.88	B.9121.36	B.9121.43
B.9120.07	B.9120.16	B.9120.31	B.9120.40	B.9120.73				

Bipolar Non-Stick US Pin

B.9129.00	B.9129.08	B.9129.18	B.9129.32	B.9129.66	B.9129.74	B.9129.82	B.9129.98	B.9229.38
B.9129.01	B.9129.10	B.9129.19	B.9129.34	B.9129.67	B.9129.76	B.9129.84	B.9229.02	B.9229.39
B.9129.02	B.9129.12	B.9129.20	B.9129.36	B.9129.68	B.9129.78	B.9129.85	B.9229.20	B.9229.40
B.9129.04	B.9129.13	B.9129.22	B.9129.37	B.9129.70	B.9129.79	B.9129.86	B.9229.26	B.9229.42
B.9129.06	B.9129.14	B.9129.30	B.9129.38	B.9129.72	B.9129.80	B.9129.88	B.9229.36	B.9229.43
B.9129.07	B.9129.16	B.9129.31	B.9129.40	B.9129.73				

Bipolar with Irrigation

B.9106.06	B.9106.14	B.9106.70	B.9106.78	B.9107.06	B.9107.12	B.9107.66	B.9107.72	B.9107.78
B.9106.08	B.9106.16	B.9106.72	B.9106.80	B.9107.08	B.9107.14	B.9107.68	B.9107.74	B.9107.80
B.9106.10	B.9106.66	B.9106.74	B.9106.82	B.9107.10	B.9107.16	B.9107.70	B.9107.76	B.9107.82
B.9106.12	B.9106.68	B.9106.76						

Monopolar

B.9105.14	B.9105.16	B.9105.19	B.9105.22	B.9105.26	B.9105.31	B.9105.36	B.9105.40	B.9105.42
B.9105.15	B.9105.18	B.9105.21	B.9105.25	B.9105.30	B.9105.33	B.9105.38		

Monopolar without plug

B.9108.14	B.9108.18	B.9108.19	B.9108.21	B.9108.22	B.9108.25	B.9108.26	B.9108.30	B.9108.31
B.9108.15								

Content

Explanation of symbols.....	4
Notice.....	5
Intended use.....	5
Contraindications.....	5
Incidents that have been reported in connection with the use of electrosurgical systems.....	5
Use and safety instructions.....	6
Reprocessing.....	6
Machine reprocessing.....	6
Manual reprocessing.....	7
Functional test and packaging.....	7
Sterilisation.....	7
Information on the validation of the reconditioning.....	8
Handling.....	8
Disposal.....	9
Guarantee and Warranty.....	9

We thank you for your confidence in our products and ask you to observe our information. Correct compliance with these instructions will ensure that you will receive optimum performance and the longest possible life for Gimmi® products.

Explanation of symbols



Caution! Refer to user manual for important safety-related information such as warnings and precautions.



Observe user manual



Manufacturer of the medical device



CE marking for European market access and Identification number of the notified body



Keep in a dry place



Order number / catalogue number



Batch designation



Non-sterile



Attention: According to US-laws, this device must only be sold by a doctor or on the instruction of a doctor.

Notice



Please read all information contained in this insert.

Incorrect handling and care, as well as misuse, can lead to premature wear of surgical instruments or risks to patients and users.

Intended use

Bipolar and monopolar forceps are used for grasping, dissection and coagulation of biological tissues. The fully assembled instrument (if assembly is needed) has to be connected – with the appropriate cable – to monopolar or bipolar output of an HF generator. Only the defined parameters have to be used.

Maximum output voltage of the generator, U_{max} :

Bipolar forceps: 500 Vp or 300 Vp

(depending on the type, see catalogue)

Monopolare forceps: 2000 Vp

Appropriate connecting cables for bipolar forceps:

Gimmi® bipolar cable / flat plug, ref. no. T.5768.XX (except T.5768.50 and T.5768.82)

Gimmi® bipolar cable / 2-pin plug, ref. no. T.5767.14 und T.5767.20

Appropriate connecting cables for monopolar forceps:

Gimmi® monopolar cable ref. no. T.9539.XX, T.9540.XX, T.9543.XX (4 mm Instrument side plug)



Instruments for electrosurgery must only be used by persons who have been specially trained or instructed in this.

Contraindications

Do not use the instrument if, in the opinion of the attending physician, the risks to the patient outweigh the benefits. The instrument is also contraindicated for cardiovascular system and central nervous system.

Incidents that have been reported in connection with the use of electrosurgical systems

- Unintended activation with resulting tissue injury in the wrong location and/or damage to the equipment.
- Fire in connection with surgical drapes and other inflammable materials.
- Alternating current paths leading to burns on spots where the patient or user comes into contact with components without insulation.
- Explosions caused by sparks in the proximity of inflammable gases.
- Perforation of organs. Sudden severe bleedings.

Use and safety instructions

Non-observance of these use and safety instructions may lead to injuries, malfunctions or other unexpected incidents.

- When using electrosurgery in patients with pacemakers or other active implants, special requirements apply (e.g. low HF-current, patient monitoring). In any case, a cardiologist or appropriate medical specialist must be consulted.
- Before initial use and any further use, all instruments must be completely cleaned, disinfected and sterilised and their function must be checked.
- It is very important to check every surgical instrument for visible damage and wear, such as cracks, breaks or insulation defects before each use. In particular areas such as blades, tips, notches, locking and blocking devices, as well as all movable parts, insulations and ceramic elements must be checked carefully.
- Never use any damaged instruments.
- Never use the instruments in the presence of flammable or explosive substances.
- When temporarily not in use, the instrument must be placed electrically insulated from the patient.
- Activate electrosurgical current only if the contact areas are in full view and have good contact with the tissue that needs to be treated. Do not touch any other metallic instruments, trocar sleeves, optics or similar objects during use.
- Observe the use and safety instructions of the manufacturer of the high-frequency surgical device.



Applies for monopolar mode of operation: Ensure correct application of the neutral electrode on the patient; otherwise, there is a danger of burns.

Reprocessing

Due to the product design, the materials used and the intended purpose, it is not possible to define a limit with regard to the maximum possible number of reprocessing cycles. The serviceable life of the instruments is determined by their function as well as by a careful handling. Instruments for electrosurgery are by their nature subject to increased wear depending on the type and time of use.

Preparation and transport

Immediately after each use, clean the instruments with a soft brush under cold tap water until all visible contamination is removed. Do not use fixation agents or hot water (>40°C). Storage and transport of the instruments to the reprocessing location must take place in a sealed container.

Machine reprocessing

Cleaning

Place the instruments in a basket on the insert module or on the inserts of the MIS module and start the cleaning process.

1. Prerinse with cold water for 1 min
2. Discharge
3. Prerinse with cold water for 3 min.
4. Discharge
5. Wash at 55°C with a 0.5% alkaline or at 45°C with an enzymatic cleaning agent for 5 min.
6. Discharge
7. Neutralise with warm tap water (>40°C) and a neutralising agent for 3 min.
8. Discharge
9. Rinse with warm tap water (>40°C) for 2 min.
10. Discharge

Disinfection

Machine-operated thermal disinfection must be carried out under observation of the national requirements regarding the Ao value (see ISO 15883).

Drying

Dry the outside of the instruments by carrying out a drying cycle of the cleaning/disinfection machine.

If necessary, manual drying may additionally be carried out using a lint-free cloth. Dry cavities by blowing with sterile compressed air.

Manual reprocessing**Ultrasonic pre-cleaning**

1. The instruments are placed in an ultrasonic bath with 0.5% enzymatic cleaning detergent and treated with ultrasound for 15 minutes at 40°C/104°F.
2. Remove the instrument and rinse them completely with cold water to remove the cleaning detergent.

Cleaning

Prepare a cleaning bath according to the manufacturer's instructions.

1. Rinse products with cold tap water (<40°C) until all visible contamination has been removed. Remove adhering dirt by using a soft brush.
2. Place products in the prepared cleaning bath so that they are completely submersed. Observe residence time according to the manufacturer's instructions.
3. Clean the instrument in the bath manually using a soft brush. Brush all surfaces several times.
4. The following step only applies to channels and the insides of tubes: Push the brush into and out of the tubes at least six times. Rinse the tubes with DI water. Repeat the procedure.
5. Rinse the products thoroughly with DI water to remove the cleaning agents without residue.

Disinfection

Prepare a disinfectant bath according to the instructions of the disinfectant manufacturer. Place the instruments in the disinfectant bath and observe the specified residence time. Rinse the products very thoroughly with DI water to remove the disinfectant without residue.

Drying

Manual drying is carried out using a lint-free cloth and sterile compressed air, in particular for drying cavities and channels.

Functional test and packaging

Perform visual inspection for cleanliness and integrity. If necessary, repeat reprocessing until the instrument is visually clean.

Packaging must comply with the ISO 11607 and EN 868 standards for packaging for sterilised instruments.

Non-stick bipolar forceps

The polished precious metal forceps tips of the Non-stick bipolar forceps may tarnish similar to silver. This does not impair function.

Bipolar forceps with irrigation

The enclosed wire insert should be always inserted in the irrigation channel, except during use and cleaning, in order to prevent clogging. The irrigation channel must be rinsed very thoroughly during cleaning. The passage has to be checked after cleaning.

Sterilisation

Sterilisation of the products with fractional pre-vacuum procedure (in accordance with ISO 13060 / ISO 17665) under observation of the respective national requirements.

- 3 pre-vacuum phases with a pressure of at least 60 mbar.
- Heating up to a sterilisation temperature of at least 132°C and at most 137°C
- Exposure time: at least 3 min.; at most 18 min.
- Drying time: at least 10 min.



If contamination with prions (CJD) is suspected, differing national guidelines are to be followed and longer holding times (i.e. 15 min.) may apply.

Storage

Sterilised instruments must be stored in a dry, clean and dust-free environment. The applicable national guidelines must be followed.

Repairs

Never attempt to perform repairs yourself. Service and repair work must only be performed by persons trained and qualified accordingly. If you have any question regarding these matters, contact either the manufacturer or your medico-technical department.



Defective products must complete the entire reprocessing process before being returned for repair.

Information on the validation of the reconditioning

The following testing instructions, materials and equipment have been used for validation:

Cleaning agents (for machine use)

- Neodisher FA by Dr. Weigert (alkaline)
- Endozime by Ruhof (enzymatic)

Cleaning agents (manual cleaning)

- Cidezime, Enzol Enzym detergent, Johnson&Johnson

Disinfectants (manual disinfection)

- Cidex OPA, Johnson&Johnson

Neutralising agent

- Neodisher Z by Dr. Weigert

Cleaning and disinfection device

- Miele Desinfector G 7735 CD
- Miele insert module E 327-06
- Miele MIS module E 450

For details, see report.

SMP GmbH # 01707011901	(machine cleaning)
MDS GmbH # 135196-10	(man. cleaning/disinfection)
Nelson Labs # 200432706-02	(sterilisation)
MDS GmbH Testbericht 084183-10	(sterilisation)

If the chemicals and machines described above are not available, the user has to validate the used process accordingly.

Handling

During transport, cleaning, care, sterilisation and storage, all surgical instruments should be handled with maximum care. This applies particularly to blades, fine tips and other sensitive areas.

Disposal

Disposal must be carried out in accordance with the respective applicable local and national laws and regulations.

Guarantee and Warranty

Gimmi® GmbH exclusively supplies tested and faultless products to its customers. All products are designed and manufactured to comply with maximum quality requirements. We refuse any liability for products which have been modified as compared to the original product, misused or handled or used improperly.

In case of repairs by companies, which are not authorised for repairs by Gimmi®, the warranty is inapplicable.

If a serious incident occurs in connection with a product of the company Gimmi®, the manufacturer must be notified immediately. After consultation with the latter, a notification is sent to the competent authority of the member state in which the user is established.

Gimmi® GmbH guarantees that this product is free of defects in material and design at the time of purchase. The warranty period is 2 year from date of purchase.

This guarantee applies to all purchases and is limited to free-of-charge repair or replacement of the product if it has been established that a material or manufacturing fault existed.

Normal wear and tear damage or misuse in the intended purposes are exempt from guarantee. Shipping costs and risk of loss or damages are not covered by Gimmi® GmbH. Prior to returning the product for maintenance purposes or for warranty claims, the product must be cleaned, disinfected, and sterilized.

Disclaimer

The manufacturer does not assume liability for immediate or secondary damages and the warranty claim is void if

- the product was improperly used, reprocessed or serviced,
- instructions and provisions stated in the operating manual were not observed,
- non-authorized personnel performed repairs, calibrations or modified the product in any manner.

The company Gimmi® assumes no liability, if this customer information is demonstrably violated.



Manufacturer

Gimmi® GmbH
Carl-Zeiss-Strasse 6
78532 Tuttlingen
Germany

contact@gimmi.de
www.gimmi.de